

	GENEESS GUIDELINE BEST PRACTICES FOR ELECTRONICS IN MEDICAL DEVICES			
0	• Version 1: May 2015			
0	Content:			
	<ul> <li>Basic elements of regulatory requirements for medical devices.</li> </ul>			
	<ul> <li>The realization of medical devices containing electronics.</li> </ul>			
	<ul> <li>Basic elements of project management for medical devices.</li> </ul>			
	• Basic elements of Design-for-X: Design-for-Test, Design-for-Reliability, Design-for-RoHS, etc.			
	<ul> <li>Basic overview of biocompatibility for medical devices.</li> </ul>			
	<ul> <li>Design &amp; testing</li> </ul>			
0	Authors:			
	<ul> <li>Maaike Op de Beeck - imec</li> </ul>			
	• Filip Ponsaerts - imec			
	<ul> <li>Frederik Horemans</li> <li>- DSP Valley</li> </ul>			
ima	C © imec 2016 / CONFIDENTIAL 7 WWW.CEDM.BE			

## **GENEESS GUIDELINE**

APPROVAL FDA (US) AND CE MARKING (EU)

- Pre-approval in US
  - FDA (Center for Devices and Radiological Health)
    - Class I: very low to low risk
    - o Class II: moderate risk
    - Class III: high risk (all implants are class III)
  - Quality System Regulation
- Approval in EU
  - Medical Devices Directive (MDD): 93/42/EEC & 2007/47/CE
    - Class I, Is, Im
    - Class IIa
    - Class IIb
    - $\circ~$  Class III: highest risk (such as implants)
  - Notified Bodies (NB) assess whether a product meets the MDD (CE marking)
  - Quality Management System
  - IEC 60601-1


© imec 2016 / CONFIDENTIAL

8

WWW.CEDM.BE

<u>Jedm</u>

<b>GENEESS GUIDELINE</b> PROJECT MANAGEMENT BEST PRACTICES		
<ul> <li>Choice of the type and level of project management</li> <li>Quality Management System (ISO 13485)</li> </ul>		
• Recommended minimum:		
<ul><li>Design Control</li><li>Define the Scope and Objectives</li></ul>		
<ul> <li>Define the Deliverables</li> <li>Project Planning</li> </ul>		
<ul> <li>Communication plan</li> </ul>		
<ul> <li>Tracking and Reporting</li> </ul>		
<ul> <li>Change Management</li> </ul>		
<ul> <li>Risk Management</li> </ul>		
imec © imec 2016 / CONFIDENTIAL 9	WWW.CEDM.BE	



ſ

<b>GENEESS GUIDELINE</b> DFX GUIDELINES	
<ul> <li>Important when designing a Printed Board Assembly</li> <li>Do's and don'ts of good DfX practice</li> <li>Basic Design-for-Assembly rules</li> <li>Basic Design-for-Manufacturing rules</li> <li>Basic Design-for-Test rules</li> <li>Basic Design-for-Reliability rules</li> </ul>	
inec © imec 2016 / CONFIDENTIAL 11	

## **GENEESS GUIDELINE**

BIOCOMPATIBILITY, BIO STABILITY AND STERILITY

- Implantable electronic devices
- Wearable electronic devices
- Natural bio-response upon implantation of foreign material
- o Biocompatibility related to material-tissue interaction
- Biocompatibility realized by dedicated encapsulation
- $\circ~$  3 most used types of sterilization:
  - Radiation based techniques
    - $\,\circ\,$  Gamma radiation
    - $\,\circ\,$  Electron beam radiation
    - X-ray radiation
  - ETO sterilization
  - Autoclave sterilization

imec 2016 / CONFIDENTIAL


12

WWW.CEDM.BE

<u>Jedm</u>


٦

GENEESS GUIDELINE ISO 13485 CHANGES MARCH 2016			
<ul> <li>Some of the changes:</li> <li>Risk-based approach beyond product realization.</li> <li>Increased linkage with regulatory requirements (regulatory documentation)</li> <li>Validation of software applications (QMS software, Process control, measurement &amp; monitoring)</li> <li>Appropriate infrastructure (sterile medical devices, clean environment,)</li> <li>Planning &amp; documenting CAPA, implement without delay</li> </ul>			
<ul> <li>Guideline will be updated accordingly</li> </ul>			
inec 2016 / CONFIDENTIAL <sup>13</sup> WWW.CEDM.BE	<u>LEDM</u>		



٦

GENEESS TOOLKIT ACTIVE MEDICAL DEVICES				
<ul> <li>Assess strengths and weaknesses</li> <li>Assessment of active medical device</li> <li>Classification based on intended use:         <ul> <li>Class I,ls,lm</li> <li>Class Ila</li> <li>Class Ilb</li> <li>Class III</li> </ul> </li> <li>First assessment of required Quality management system</li> <li>Links, documents &amp; checklist</li> <li>Evaluate impact of intended use :         <ul> <li>On classification</li> <li>On CE marking</li> </ul> </li> </ul>	<ul> <li>Assessment of active medical device</li> <li>Classification based on intended use: <ul> <li>Class I,Is,Im</li> <li>Class IIa</li> <li>Class IIb</li> <li>Class III</li> </ul> </li> <li>First assessment of required Quality management system</li> <li>Links, documents &amp; checklist</li> <li>Evaluate impact of intended use : <ul> <li>On classification</li> </ul> </li> </ul>			
imec 2016 / CONFIDENTIAL 15				



GENEESS TOOLKIT ACTIVE MEDICAL DEVICES			
<ul> <li>Intended use: measure heart rate</li> </ul>			
<ul> <li>Diagnostic use e.g. 24 hour follow up</li> </ul>			
<ul> <li>Follow up of a heart condition as part of or after treatment</li> </ul>			
<ul> <li>Clinical trial</li> </ul>			
• Stress test:			
<ul> <li>Diagnostic</li> </ul>			
• Prevention			
○ Study			
<ul> <li>Sport environment:</li> </ul>			
<ul> <li>Early detection of risk factors</li> </ul>			
<ul> <li>Training progress</li> </ul>			
inec 2016 / CONFIDENTIAL 17	WWW.CEDM.BE	<u>LEDM</u>	

<b>GENEESS T</b> WEB-BASED	OOLKIT		
	Geneess Toolkit voo x	(A) - 0 - 3	4
	← → C D localhost/Geneesstoolkit/toolkit.php?toolkitnr=0	\$☆© ≣	4
	Actieve medische hulpn Toolkit voor KMO's en		
		EN & Waanderen	
	De todiki Met deze toolkit beogen wy de KMO's en Skatten alvast op de weg te helpen om hun meiskede en of zaspesterste tote ty Deze toolkit ped de mogelijkheid om een erete inschaften of de hundlichte neu of zaspesterste tote ty Deze toolkit ped de mogelijkheid om een erete inschaften de hundlichte neu opzahle, een erete De konte behandel alten en de toolkit behandel alten en ander het en genaphietende streve medicke hund De konte behandel alten de stand genaphiet om hundlichte interver medicke hundlichte interver medicke hund De konte behandel alten de stand genaphiet om hundlich hundlichte interver medicke hund De konte behandel alten de hundlicht hundlichte hundlich hundlich hundlichte medick hundlich hundlichte medick hundlichte hundlichte hundlichte hundlichte hundlichte hundlichte medicken hundlichte hundlitte hundlichte hundlitte hundlichte hundlichte hu	assen. isch halpmiddel op de Europese markt te zetten. (CE markering van MDD) teolichtung van de reglementering en handige hjutjes. middelen exclusief in-vitro halpmiddelen	
	<u>GENEESS</u> project mast <u>cEDMinus</u> ; en [ GEzonder door Nano-Elektronica En Slimm	SP Valley 1 Specialisatie	
			0
imec © imec 2016 / CONF	IDENTIAL 18		WWW.CEDM.BE

